
www.lagis.com.tw/eifu_en.php

BEFORE USING PRODUCT, READ THE FOLLOWING INFORMATION THOROUGHLY.

◆ DEVICE DESCRIPTION:

The LAGIS Balloon Trocar is a single-use device sterilized by Gamma Radiation. It is composed of a cannula with balloon, an obturator, a stopcock and a syringe. The diameter of the cannula is 12 mm, which is able to accommodate endoscopic instruments from the size of 5mm to 12mm. In order to be used in different patients and surgery, the length of the device comes in two sizes: Standard and Extra-long. The fixed components of the cannula also come in two types: Hasson Cone and Stability Ring. The obturator with optical tip has a clear and transparent tip. It is designed to offer direct visualization with the endoscope during insertion. The cannula contains a set of seal to minimize the risk of air leakage when endoscopic instruments are inserted or withdrawn through the LAGIS Balloon Trocar. A stopcock is located at the side of the cannula for the connection to the luer lock to provide gas insufflation and desufflation.

◆ INDICATIONS/INTENDED PURPOSE:

- The LAGIS Balloon Trocar can be applied to general, thoracic, gynecological, laparoscopic, abdominal and other minimally invasive surgical procedures to penetrate the body cavity and establish a path of entry for endoscopic instruments.
- Expected clinical benefits: reduce the surgery risks, shorten the surgery time, minimize the wound, and reduce the bleeding situation.
- Patient target group: the LAGIS Balloon Trocar is intended for adults who need minimally invasive procedures.

◆ CONTRAINDICATIONS:

The LAGIS Balloon Trocar is not intended for use when minimally invasive techniques are contraindicated.

◆ INSTRUCTIONS FOR USE:


Prior to use, inspect the contents of the packaging and the packaging itself to ensure the integrity is not compromised.

- Verify the balloon is fully deflated prior to use. To deflate the balloon, attach the syringe onto the balloon inflation valve and pull out the syringe.
- Insert the obturator into the cannula until the tip is exposed and make sure it is fully assembled.
- Make an incision at the puncture site just long enough to accept the LAGIS Balloon Trocar to penetrate the body cavity.
- Use the provided syringe to inflate the balloon through the balloon inflation valve with 10cc of air.
- Slide the Hasson Cone down into the incision. Lock the Hasson Cone into place by pressing the button in the Locked position. If desired, pull and wrap the sutures around the suture slots to secure the Hasson Cone. (TRB-12V1, TRB-12VX1); Slide the Stability Ring on the incision (TRS-12V1, TRS-12VX1).
- Retract the obturator out of the cannula once the LAGIS Balloon Trocar is securely installed.
- Perform the surgery according to the endoscopic surgical procedures.
- To remove the cannula upon completion of the procedure, attach the syringe onto the balloon inflation valve and pull out the syringe to deflate the balloon. Unwind sutures from the Hasson Cone. Remove the LAGIS Balloon Trocar.

9. Dispose of the used LAGIS Balloon Trocar properly and safely by healthcare professionals in accordance with the national waste management regulation.

◆ WARNINGS AND PRECAUTIONS:

- Do not use if the device or the sterile packaging is open or damaged, as this may cause an increase in the incidence of wound infections.
- The LAGIS Balloon Trocar is designed to be used with rigid endoscopes. Do not use the LAGIS Balloon Trocar with flexible endoscopes.
- If endoscopic instruments or accessories from different manufacturers are used together in a procedure. Please verify compatibility prior to performing the procedure.
- Using an endoscopic instrument with a diameter which is smaller than the specification specified for the LAGIS Balloon Trocar may result in desufflation of the body cavity.
- The surgical site should be inspected for hemostasis after removing the LAGIS Balloon Trocar from the body cavity.
- This device is packaged and sterilized for single patient use only. Do not reuse, reprocess, or resterilize. Reuse, reprocess, or resterilization may create contamination, infection, cross-infection, or even compromise the integrity of the device and lead to device failure which in turn may result in patient injury, illness or death.
- This device should be operated only by surgeons with adequate training and familiarity with surgical techniques. Consult medical literature for further information regarding techniques, hazards, contraindications, and complications prior to performing the procedures.
- Return the device and packaging to the distributor or implement safe disposal handling by the healthcare professionals if the sterile packaging is damaged or unintentionally opened before use.
- The used device which comes into contact with bodily fluids should be considered as biomedical waste and safely disposed of by healthcare professionals in compliance with the national waste management regulation to prevent patients, users, and other persons from biological contamination or infection.
- In the event of malfunction of the device, changes in performance that they may affect safety, or even any serious incident. Please report to the manufacturer, the competent authority of the Member State and the authorized representative immediately according to the contact information provided on the label or in the instructions for use.
- If the user feels the resistance is higher than usual when inserting the endoscopic instruments, the angle of insertion should be adjusted to align with the axis of the LAGIS Balloon Trocar for fear of damaging the valve. (Please see the figure on the right-hand side.)
- Improper assembly or operation of the device may result in tissue or organ injuries from mechanical damage to the patient or operator.
- Excessive inflation of the balloon can cause balloon to burst.
- Do not allow sharp instruments or objects to come into contact with the balloon. Please be careful when using the device around metal clamps, sutures, and staples.
- Ensure insufflation tubing is connected to insufflation port on Stopcock.
- Inflating the balloon can cause rupture when the balloon does not completely reach the abdominal cavity.

◆ ENVIRONMENTAL CONDITIONS FOR STORAGE:

Appropriate storage environment is a clean and dry area away from sunlight with a temperature range of 13~30°C (55.4~86°F).

◆ EXPLANATION OF SYMBOLS :

	Temperature limit		Caution		Authorized representative in the European Community/European Union
	Do not resterilize		Manufacturer		Not made with natural rubber latex
	Do not re-use		Consult instructions for use or consult electronic instructions for use		Do not use if package is damaged and consult instructions for use
	Date of manufacture		Single sterile barrier system		Sterilized using irradiation
	Medical device		Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.		
	CE Marking according to MDR 2017/745		UKCA Marking according to UK MDR 2002		

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